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The electronic baby scale has been designed to help you monitor the growth of your baby, accurately and reliably quantifying how much milk has been assimilated at each feed. The TARE function allows you to weigh the baby using a soft cloth or towel.

The scale come with a useful WEIGHT-BLOCK function that allows a steady weight to be displayed despite the child wriggling about. Taking up to 20 kg, and with 10g divisions, the scales will also help you monitor your child's progress in terms of weight gain during all stages of growth. As soon as the child can stand up on his own, weighing operations can be performed without the pan, transforming the unit into practical bathroom scale for your infant. The new height-measurement function allows you to monitor the growth in terms of height of your baby.

Laica S.p.A. - Viale del Lavoro, 10 - Fraz. Ponte - 36021 Barbarano Vicentino (VI) - ITALY  
Fax: +39-0444-795324

Made in China

€ 0434 (Dir. 93/42/EEC)

**YOU ARE STRONGLY RECOMMENDED TO READ THE INSTRUCTIONS AND WARNINGS HEREIN BEFORE USING THE EQUIPMENT AND TO KEEP THE BOOKLET IN A SAFE PLACE.**

#### TIPS FOR USE

- **Place the scale on a flat and hard surface.**  
A soft and uneven floor (e.g. carpets, linoleum) could cause an alteration to the weight reading.
- **If the scale is used without plate, the baby has to stand on the scale with the feet parallelly positioned and the weight must be equally distributed (see the drawing), before the zeros disappear. Pay attention the baby stands still while the scale measures the weight. Pay attention that the baby does not lean on anything while being measured, as this could alter the reading.**
- **During the measurement of height, the child must be kept monitored by an adult.**

#### INSTALLING BATTERIES

Use a 9V alkaline battery only. Remove the battery cover on the underside of the scale. Insert the battery and refit the cover. Remove the run down batteries by lifting the cover of the battery compartment located on the bottom of the device. Extract the batteries and dispose them as special waste at the specific collection points forecast for recycling. For more information about the disposal of run down batteries, please contact either the store where the device containing the batteries was purchased, or the Municipality, or the waste disposal local service. **Caution: For greater safety, children younger than 12 years should not remove the batteries.**

#### FUNCTION KEYS

##### ON

The ON key turns the scale on.

##### OFF

The OFF key turns the scale off.

##### SELECT UNIT (UNIT)

Press this key to select the unit of measure you require: kg or pounds. These scale has a capacity of 20 kg, division 10 g. or 44.40 pounds, division 1/2 oz.

After turning the scale on, if the display shows 0.00 the weight is expressed in kg and g, if the display shows 0:0 the weight is expressed in pounds and ounces.

##### HOLD (BLOCK WEIGHT FUNCTION)

Due to the baby moving, the weight shown on the display changes all the time. When the baby is laid on the plate, press the "HOLD" key. The text "hold" blinking on the display indicates that the scale is calculating automatically the average of the detected weights.

When the text "hold" stops blinking on the display, the weight of the child is displayed.

##### TARE

If the baby is laid on a towel or cover, use the TARE function.

Turn the scale on with the ON key.

The display shows "8888". Wait until the display shows "0.00" then place the towel or cover on the plate. Press the TARE key and wait until the display shows "0.00"; the scale is now ready for use.

##### MODE

The scale is equipped with the new height-measurement function.

When the child can stand safely alone, it is possible to measure his/her height to monitor his/her growth.

Remove the plate from the scale. Switch on the scale by pressing the key ON. The display shows "8888".

Wait until the display shows "0.00". Let the child get on the platform as described above and press the key MODE. On the display appears "0.0 cm".

Pull upwards for 2/3 cm the knob positioned over the display and insert it in the height fixing accessory according to figure A.

Pull the accessory upwards slowly and without sudden movements until it is positioned on the head of the child (figure B).

Now, the display shows the height in cm; the text "hold" indicates that the data is definitive (average of the perceived heights). It is possible to measure the height in feet by pressing the key UNIT before measuring the height of the child or just after the data has appeared on the display.

After having ended the height measurement procedure, accompany slowly the knob.

#### USING THE SCALE WITH THE PLATE

Remove the 4 red rubber plugs.

Fix the plate to the scale using the supplied locking nut (see the illustration).

## Instructions and guarantee

Press the ON key.

The display shows "8888". Wait until the display shows "0.00"; the scale is now ready for use. Lay the baby on the plate and wait until the display shows its weight.

Press the OFF key to turn the scale off, otherwise the scale turn off automatically after 2 minutes without being used.

#### USING THE SCALE WITHOUT THE PLATE

This baby scale can also be used without the plate, when the baby can stand up on his own and does not weigh more than 20 kg.

Remove the plate and keep the locking nut for future use.

Replace the 4 red rubber plugs. Turn the scale on and wait until the display shows "0.00".

Get the baby to stand on the scale.

#### MAINTENANCE

- Take care of your scale, as it is a precision instrument. Therefore, avoid exposing to extreme temperatures, damp, knocks, dust and direct sunlight. Do not drop or knock the scale.
- Keep the scale away from sources of heat.
- **Do not jump on the scale**
- Clean the platform with a damp cloth
- Clean the bottom of the scale with just a soft dry cloth.
- Never use chemical or abrasive products.
- Take care that no liquids penetrate into the scale, as it is not waterproof.
- Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquids could cause corrosion and effect the precision of the scale.

#### SYMBOLS

- "LO" indicates that the batteries are flat. Replace them.
- "ERR" shows the weight is too much.
- The scale is ready for use when the display shows "0.00"; if the instructions have not been followed correctly, the display shows "ERR" or "----".

#### TECHNICAL DATA

Range of measurement: 1 kg - 20 kg / 30 cm - 150 cm

Division: 10 g

Battery: 1 alkaline battery 6LR61 series

Accuracy within: +/-1% + 1 division of applied weight

Working conditions: +10°C +40°C RH 85%

Storage conditions: +10°C 50°C RH 85%

Attention! Read carefully the instructions for use.

#### DISPOSAL PROCEDURE (DIR. 2002/96/EC - WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

#### WARRANTY

This equipment is covered by a 2-year warranty elapsing from the purchase date, which has to be proven by means **of the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force (Italian legislative decree No. 24 dated Feb. 02, 2002 and following "Consumers' Act" Italian legislative decree No. 206 dated Sep. 03, 2005) and applies only in case the consumer is a private entity. Laica products are designed for domestic use; therefore, their use in public businesses is forbidden.

The warranty covers only production defects and is void, if the damage is caused by accidental events, wrong use, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories may lead to warranty voidness.

Do not open the equipment for any reason whatsoever; in case of opening or tampering of the equipment, warranty decays definitively. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are standard supplied. Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com. Nothing shall be owed for the repairing or replacement of products covered by the warranty terms. In case of faults, refer to a dealer; DO NOT send anything directly to LAICA. All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.

The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the warnings related to installation, use and maintenance of the equipment. Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

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